



Anti-Mouse IgG1 polyclonal antibody [HRP] (CPBT-68024GM)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Product Overview	Goat anti Mouse IgG1 antibody recognizes Mouse IgG1. This antibody has been cross absorbed against mouse IgM, IgG2a, IgG2b, IgG3 and IgA, pooled human sera and purified human paraproteins. Goat anti Mouse IgG1 antibody shows minimal cross-reactivity with human immunoglobulins.
Specificity	IGHG1
Immunogen	Mouse IgG1 paraproteins.
Isotype	IgG
Source/Host	Goat
Species Reactivity	Mouse
Conjugate	HRP
Applications	IHC-Fr; ELISA; IHC-P; WB
Format	Purified IgG conjugated to Horseradish Peroxidase (HRP) - liquid
Size	1 ml
Buffer	50% Glycerol/50% Phosphate buffered saline, pH 7.4
Preservative	None
Storage	in frost-free freezers is not recommended. Avoid repeated freezing and thawing as this may denature the antibody. Should this product contain a precipitate we recommend microcentrifugation before use.

GENE INFORMATION

Gene Name	Ighg1 immunoglobulin heavy constant gamma 1 (G1m marker) [Mus musculus (house mouse)]
Official Symbol	IGHG1
Synonyms	IGHG1; immunoglobulin heavy constant gamma 1 (G1m marker); IgG1; Igh-4; VH7183;
Entrez Gene ID	16017
UniProt ID	P01868
Chromosome Location	12 F1; 12 62.1 cM
Function	antigen binding;