



Anti-Mouse IgM polyclonal antibody [R-PE] (CPBT-68047GM)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Product Overview	Goat anti Mouse IgM antibody recognizes murine IgM and is cross-adsorbed to reduce reactivity with other mouse heavy or light chains. Flow Cytometry Use 50ul of the suggested working dilution to label 106 cells in 100ul.
Specificity	IGHM
Immunogen	Mouse IgM mu heavy chain.
Isotype	IgG
Source/Host	Goat
Species Reactivity	Mouse
Conjugate	PE
Applications	FC
Format	Purified IgG conjugated to R. Phycoerythrin (RPE) - lyophilised
Size	500 µl
Preservative	0.09% Sodium Azide
Storage	Prior to reconstitution store at +4°C. Following reconstitution store at +4°C. This product should be stored undiluted. DO NOT FREEZE. This product is photosensitive and should be protected from light. Should this product contain a precipitate we recommend microcentrifugation before use.

GENE INFORMATION

Gene Name	Ighm immunoglobulin heavy constant mu [Mus musculus (house mouse)]
Official Symbol	Ighm
Synonyms	IGHM; immunoglobulin heavy constant mu; Igm; muH; Igh6; Igh-6; Igh-M; AI326478; Ig mu; immunoglobulin heavy chain 6 (heavy chain of IgM); mu-immunoglobulin;
Entrez Gene ID	16019
UniProt ID	P01873
Chromosome Location	12 62.1 cM; 12 F1-2
Pathway	Adaptive Immune System; Antigen activates B Cell Receptor (BCR) leading to generation of second messengers; Immune System; Signaling by the B Cell Receptor (BCR);
Function	antigen binding; protein binding; transmembrane signaling receptor activity