



Rabbit Anti-Ganglioside GD1b polyclonal antibody (CPBT-LL010)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Specificity	No cross-reaction with other carbohydrate epitopes
Target	Ganglioside GD1b
Immunogen	Purified ganglioside GD1b and complete Freund's adjuvant was used to immunize rabbits
Isotype	IgG/IgM
Source/Host	Rabbit
Species Reactivity	N/A
Purification	Serum, Unpurified
Conjugate	Unconjugated
Applications	ELISA,TLC
Molecular Weight	N/A
Format	Liquid
Size	50 µl
Buffer	Dilution: Phosphate buffered saline (pH 7.4) is recommended.
Preservative	None
Storage	-20°C
Ship	Dry ice