



Human Anti-Human Parathyroid hormone-related protein monoclonal antibody, clone 2F8-10-3 (CPBT-LL019)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Product Overview	Recombinant monoclonal antibody to Parathyroid hormone-related protein.
Specificity	This antibody binds human parathyroid hormone. Parathyroid hormone is a hormone secreted by the parathyroid glands that regulates the serum calcium concentration through its effects on bone, kidney, and intestine.
Target	PTHrP
Immunogen	The original antibody was generated by immunizing a human antibody-producing transgenic mouse with human PTHrP.
Isotype	IgG2
Source/Host	Human
Species Reactivity	Human
Clone	2F8-10-3
Purification	Protein A affinity purified.
Conjugate	Unconjugated
Applications	ELISA
Format	Liquid
Concentration	Lot specific

Size	200 µg, 1 mg
Buffer	PBS with 0.02% Proclin 300.
Preservative	View individual datasheet
Storage	Store at 4°C for up to 3 months. For longer storage, aliquot and store at - 20°C.
Ship	Wet ice

GENE INFORMATION

UniProt ID	P12272
-------------------	------------------------