



# Human Anti-Human sFlt-1 monoclonal antibody, clone HFB1 (CPBT-LL022)

This product is for research use only and is not intended for diagnostic use.

## PRODUCT INFORMATION

<b>Specificity</b>	sFlt-1
<b>Target</b>	sFlt-1
<b>Source/Host</b>	Human
<b>Species Reactivity</b>	Human
<b>Clone</b>	HFB1
<b>Purification</b>	Affinity purified
<b>Conjugate</b>	Unconjugated
<b>Applications</b>	ELISA(Cap) We recommend the following for sandwich ELISA (Capture - Detection): CPBT-LL022 - CPBT-LL023
<b>Format</b>	Liquid
<b>Concentration</b>	Lot specific
<b>Buffer</b>	PBS
<b>Preservative</b>	0.02% Sodium Azide
<b>Storage</b>	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze-thaw cycles.
<b>Ship</b>	Wet ice