



Human Anti-Human sFlt-1 monoclonal antibody, clone HFB2 (CPBT-LL023)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Specificity	sFlt-1
Target	sFlt-1
Source/Host	Human
Species Reactivity	Human
Clone	HFB2
Purification	Affinity purified
Conjugate	Unconjugated
Applications	ELISA(Det) We recommend the following for sandwich ELISA (Capture - Detection): CPBT-LL022 - CPBT-LL023
Format	Liquid
Concentration	Lot specific
Buffer	PBS
Preservative	0.02% Sodium Azide
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze-thaw cycles.
Ship	Wet ice