



# Human Anti-Human PLGF monoclonal antibody, clone HPB1 (CPBT-LL024)

This product is for research use only and is not intended for diagnostic use.

## PRODUCT INFORMATION

|                           |                                                                                                               |
|---------------------------|---------------------------------------------------------------------------------------------------------------|
| <b>Specificity</b>        | PLGF                                                                                                          |
| <b>Target</b>             | PLGF                                                                                                          |
| <b>Source/Host</b>        | Human                                                                                                         |
| <b>Species Reactivity</b> | Human                                                                                                         |
| <b>Clone</b>              | HPB1                                                                                                          |
| <b>Purification</b>       | Affinity purified                                                                                             |
| <b>Conjugate</b>          | Unconjugated                                                                                                  |
| <b>Applications</b>       | ELISA(Cap)<br>We recommend the following for sandwich ELISA (Capture - Detection):<br>CPBT-LL024 - CPBT-LL025 |
| <b>Format</b>             | Liquid                                                                                                        |
| <b>Concentration</b>      | Lot specific                                                                                                  |
| <b>Buffer</b>             | PBS                                                                                                           |
| <b>Preservative</b>       | 0.02% Sodium Azide                                                                                            |
| <b>Storage</b>            | Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze-thaw cycles.                      |
| <b>Ship</b>               | Wet ice                                                                                                       |