



3-(4-Hydroxy-3-nitrophenyl acetamido) propionic acid [HRP] (DAG-15365R)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Product Overview	3-(4-Hydroxy-3-nitrophenyl acetamido) propionic acid conjugated with HRP
Specificity	This product is ELISA-paired with DMAB15365
Species	N/A
Conjugate	HRP
Format	Liquid
Concentration	Batch dependent - please inquire should you have specific requirements.
Size	1 mg
Buffer	0.01M PBS, pH7.4
Preservative	None
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid freeze / thaw cycles.