



HCV non structural protein 3 [His] (DAG-2800)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Antigen Description	HCV NS3 is the first HCV antigen, which has been used in HCV diagnostic kits.
Species	HCV
Purity	>95% pure, SDS Gel Electrophoresis Result
Conjugate	His
Molecular Weight	NS3~33kD
Format	Purified, Liquid
Concentration	7.0 mg/ml
Size	1 mg
Buffer	PBS, pH 7.4±0.1
Preservative	None
Storage	Store at 4 °C. Avoid multiple freeze/thaw cycles.

BACKGROUND

Keywords	HCV NS3 transactivated protein; NS 3; NS3; NS3P; p70; Serine protease/NTPase/helicase; Hepatitis C Virus NS3; Flaviviridae; Hepacivirus; Hepatitis C virus; HCV NS-3; HCV NS3 Genotype 1a; Hepatitis C Virus NS3 Genotype 1a
----------	--