



HCV non structural protein 5 [His, Trx] (DAG-2801)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Antigen Description	HCV NS5 is an important HCV antigen.
Species	HCV
Purity	>95% pure, SDS Gel Electrophoresis Result
Conjugate	His, Trx
Molecular Weight	NS5~50kD
Format	Purified, Liquid
Concentration	5.0 mg/ml
Buffer	PBS, pH 7.4±0.1
Preservative	None
Storage	Store at 4 °C. Avoid multiple freeze/thaw cycles.

BACKGROUND

Keywords	HCV NS5; Hepatitis C virus non structural protein 5; Hepatitis C virus nonstructural protein 5; Hepatitis C Virus NS5; HCV NS5 Genotype 2a; Hepatitis C Virus NS5 Genotype 2a
----------	---