



## HCV Core [His] (DAG-2802)

This product is for research use only and is not intended for diagnostic use.

### PRODUCT INFORMATION

Antigen Description	HCV NS3 has been used in HCV diagnostic kits.
Species	HCV
Purity	>95% pure, (SDS PAGE Non-reducing)
Conjugate	His
Molecular Weight	Core~15kD
Format	Purified, Liquid
Concentration	5.0 mg/ml
Buffer	PBS, pH 7.4±0.1
Preservative	None
Storage	Store at 4 °C. Avoid multiple freeze/thaw cycles.

### BACKGROUND

Keywords	HCcAg; Core protein p19; HCV core antigen; HCV core protein; Hepatitis C Virus core protein; HCV-1 Core Ag; Hepatitis C Virus Core Antigen, genotype 6a; Flaviviridae; Hepacivirus
----------	--