



# Cytomegalovirus (DAG-2806)

This product is for research use only and is not intended for diagnostic use.

## PRODUCT INFORMATION

Species	CMV
Purity	Nickel column
Conjugate	Unconjugated
Format	Liquid
Concentration	2mg/ml
Size	1 mg
Buffer	0.5M NaCl, 0.02M Tris-HCl
Preservative	None
Storage	Store at -20 °C. Avoid multiple freeze/thaw cycles and light.

## BACKGROUND

Introduction	Cytomegalovirus is a member of the herpes virus group, which includes herpes simplex virus types 1 and 2, varicella zoster virus (which causes chicken pox), and Epstein Barr virus (which causes infectious mononucleosis). These viruses share a characteristic ability to remain dormant within the body over a long period.
Keywords	55 kDa immediate-early protein 1; 45 kDa immediate-early protein 2; CMV; Envelope glycoprotein B; gB; Glycoprotein B; HHV 5; HHV5; Human herpesvirus 5; IE1; IE2; UL122; UL122 regulatory protein IE2; UL123; UL123 regulatory protein IE1; UL55; UL55 envelope glycoprotein B; Cytomegalovirus; CMV