



# Recombinant *Hypoderma lineatum* Hypodermin-A (DAG-BY02)

This product is for research use only and is not intended for diagnostic use.

## PRODUCT INFORMATION

<b>Species</b>	Hypoderma lineatum
<b>Purity</b>	> 85% as determined by SDS-PAGE
<b>Conjugate</b>	His
<b>Applications</b>	Immunoassays
<b>Molecular Weight</b>	28.8 kDa
<b>Format</b>	Liquid or Lyophilized powder
<b>Concentration</b>	Batch dependent - please inquire should you have specific requirements.
<b>Size</b>	20 µg, 100 µg, 1 mg
<b>Buffer</b>	Tris-based buffer, 50% glycerol
<b>Storage</b>	Store at -20°C to -80°C

## BACKGROUND

<b>Keywords</b>	Hypoderma lineatum; Hypodermin-A; HA
-----------------	--------------------------------------

## GENE INFORMATION

<b>UniProt ID</b>	P35587
-------------------	--------