



Creatinine-HRP (DAG-CL005)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Product Overview	Creatinine Conjugate for use in immunoassays
Conjugate	HRP
Applications	ELISA, Immunoassays
Size	0.5 ml
Buffer	Supplied as a concentrate. Dilute as required and use working strength conjugate immediately after dilution
Preservative	None
Storage	Store at 2-8 °C for up to 3 months or -20 °C for long term storage

BACKGROUND

Introduction	Creatinine is produced as the result of normal muscle metabolism.
Keywords	Creatinine-HRP Conjugate, Creatinine Conjugate, Creatinine HRP, Creatinine HRP Conjugate