



Helicobacter pylori protein (DAG-CL008)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Antigen Description	Helicobacter pylori is a Gram-negative, microaerophilic bacterium that can inhabit various areas of the stomach, particularly the antrum. It causes a chronic low-level inflammation of the stomach lining and is strongly linked to the development of duodenal and gastric ulcers and stomach cancer.
Purity	Highly pure
Conjugate	Unconjugated
Format	Liquid
Concentration	Batch dependent - please inquire should you have specific requirements.
Size	1 mg
Buffer	PBS buffer
Preservative	None
Storage	Store at -20°C or lower. Avoid repeated freeze/thaw cycles

BACKGROUND

Keywords H pylori protein