



# S. scrofa Legumain, LGMN, AEP [His] (DAG-H10354)

This product is for research use only and is not intended for diagnostic use.

## PRODUCT INFORMATION

|              |  |
|--------------|--|
| Species      | S. scrofa  |
| Purity       | > 95 % as determined by SDS-PAGE   |
| Conjugate    | His  |
| Size         | 10 µg, 20 µg   |
| Preservative | None   |
| Storage      | Store it under sterile conditions at -70 °C. It is recommended that the protein be aliquoted for optimal storage. Avoid repeated freeze-thaw cycles. |

## BACKGROUND

|              |  |
|--------------|--|
| Introduction | No Information.  |
| Keywords     | LGMN; legumain; protease, cysteine, 1 (legumain) , PRSC1; LGMN1; cysteine protease 1; protease, cysteine 1; asparaginyl endopeptidase; protease, cysteine, 1 (legumain); AEP; PRSC1; |