



Recombinant HAdV-E fiber [His] (DAG-H10636)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Product Overview	A DNA sequence encoding the Human mastadenovirus E (HAdV-E) fiber (Met1-Gln425) was expressed with a polyhistidine tag at the C-terminus
Purity	> 90% as determined by SDS-PAGE
Conjugate	His
Applications	Immunoassays
Molecular Weight	46.56 kDa
Format	Lyophilized
Size	100 µg
Buffer	Lyophilized from sterile 20mM PBS, pH 7.5, 300mM NaCl, 10% glycerol
Storage	Store at -20°C to -80°C
Ship	Wet ice

BACKGROUND

Introduction	Mastadenovirus is a genus of viruses in the family Adenoviridae. Humans and other mammals serve as natural hosts. There are 51 species in this genus. The genus as a whole includes many very common causes of human infection, estimated to be responsible for 2 to 5% of all respiratory infections, as well as gastrointestinal and eye infections. Symptoms are usually mild.
Keywords	HAdV; Human mastadenovirus; HAdV-E fiber

GENE INFORMATION

