



Recombinant HAdV-B encapsidation protein IVa2 [His] (DAG-H10637)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Product Overview	A DNA sequence encoding the Human mastadenovirus B (HAdV-B) encapsidation protein IVa2 (Thr3-Lys448) was expressed with a polyhistidine tag at the N-terminus
Purity	> 95% as determined by SDS-PAGE
Conjugate	His
Applications	Immunoassays
Molecular Weight	52.06 kDa
Format	Lyophilized
Size	100 µg
Buffer	Lyophilized from sterile 20mM PBS, pH 7.0, 300mM NaCl, 10% glycerol
Storage	Store at -20°C to -80°C
Ship	Wet ice

BACKGROUND

Introduction	Mastadenovirus is a genus of viruses in the family Adenoviridae. Humans and other mammals serve as natural hosts. There are 51 species in this genus. The genus as a whole includes many very common causes of human infection, estimated to be responsible for 2 to 5% of all respiratory infections, as well as gastrointestinal and eye infections. Symptoms are usually mild.
Keywords	HAdV; Human mastadenovirus; HAdV-B IVa2

GENE INFORMATION

UniProt ID

YP_002213830.1
