



Mouse IgG [PE/TXRD] (DAG-IC109)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Isotype	IgG
Source/Host	Mouse
Species Reactivity	N/A
Conjugate	PE/TXRD
Applications	FC, ELISA, FLISA, IHC, IP, Blocking, Control
Format	Liquid
Concentration	Batch dependent - please inquire should you have specific requirements.
Size	100 µg
Buffer	Supplied in PBS/NaN3.
Preservative	See individual product datasheet
Storage	Store at 2-8°C

BACKGROUND

Introduction

Isotype controls are primary antibodies that lack specificity to the target, but match the class and type of the primary antibody used in the application. Isotype controls are used as negative controls to help differentiate non-specific background signal from specific antibody signal. Depending upon the isotype of the primary antibody used for detection and the target cell types involved, background signal may be a significant issue in various experiments.