



Hyaluronic acid [KLH] (DAG-IV05)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Specificity	N/A
Species	N/A
Conjugate	KLH
Applications	User optimized
Format	Liquid
Concentration	Batch dependent - please inquire should you have specific requirements
Size	1 mg
Preservative	None
Storage	Long Term: -20°C. Avoid repeated freezing and thawing.
Ship	Wet ice

BACKGROUND

Introduction	Hyaluronic acid is a mucopolysaccharide which is an integral part of the gel-like substance of animal connective tissue that serves as a lubricating and protective agent.
Keywords	hyaluronic acid from bovine vitreous humor; Hyaluronic acid from Cockscomb; hyaluronic acid sodium from*streptococcus zooepid; Hyaluronic Acid Food /Cosmetic Grade /Injection 90.0% Min; HYALURONICACID(HA)38%(BULK; Hyaluronic acid