



Recombinant Human PIIINP (DAG-IV06)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Species	Human
Purity	> 90 % as determined by SDS-PAGE.
Applications	ELISA, CLIA
Molecular Weight	64 kDa
Format	Liquid
Size	1 mg
Buffer	250mM Glycine, 25mM Tris-HCl, 1/1000 SDS
Preservative	None
Storage	Store at -20°C.

BACKGROUND

Introduction	PIIINP is the amino terminal peptide of type III procollagen, released during the synthesis and deposition of type III collagen. PIIINP in the serum can be derived from the synthesis of new type III collagen or from the degradation of existing type III collagen fibrils.
Keywords	PIIINP; P3NP; P III NP; N-Propeptide Of Type III Procollagen; Procollagen III Amino Terminal Propeptide