



Fumonisin [BSA] (DAG-L3200MB)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Antigen Description	The fumonisins are a group of mycotoxins derived from Fusarium, Liseola section. They have strong structural similarity to sphinganine, the backbone precursor of sphingolipids.
----------------------------	--

Species	None
----------------	------

Purity	Dialysis followed by column purification
---------------	--

Conjugate	BSA
------------------	-----

Applications	Verification by ELISA
---------------------	-----------------------

Format	Liquid
---------------	--------

Concentration	Batch dependent - please inquire should you have specific requirements.
----------------------	---

Size	1 mg
-------------	------

Buffer	0.01M pH7.4 PBS
---------------	-----------------

Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid freeze / thaw cycles.
----------------	---