



Methadone Metabolite Antigen [BSA] (DAG-NS001)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Species	N/A
Conjugate	BSA
Applications	ELISA
Format	Liquid
Buffer	PBS or Tris-Gly
Preservative	None
Storage	Store at 4°C short term. For long term storage, store at -20°C, avoiding freeze/thaw cycles.

BACKGROUND

Introduction	2-Ethylidene-1,5-dimethyl-3,3-diphenylpyrrolidine (EDDP) is a major metabolite of methadone.
Keywords	2-Ethylidene-1,5-dimethyl-3,3-diphenylpyrrolidine; EDDP; Methadone metabolite