



Dephenylethyl fentanyl Antigen [BSA] (DAG-NS002)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Species	N/A
Conjugate	BSA
Applications	ELISA
Format	Liquid
Buffer	PBS or Tris-Gly
Preservative	None
Storage	Store at 4°C short term. For long term storage, store at -20°C, avoiding freeze/thaw cycles.

BACKGROUND

Keywords	Dephenylethyl fentanyl; NFYL
----------	------------------------------