



CAPN11 peptide (DAG-P0244)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Antigen Description	Calpains constitute a family of intracellular calcium-dependent cysteine proteases. There are eight members in this superfamily. They consist of a variable 80 kDa subunit and an invariant 30 kDa subunit. This calpain protein appears to have protease activity and calcium-binding ability. A similar mouse protein may play a functional role in spermatogenesis and in the regulation of calcium-dependent signal transduction events during meiosis. [provided by RefSeq, Dec 2008]
Purity	> 95 % by SDS-PAGE.
Conjugate	Unconjugated
Applications	ELISA, WB
Format	Liquid
Buffer	Preservative: None Constituents: 0.001% Tween 20, 30mM HEPES, 2mM EDTA, 150mM Sodium chloride, pH 6.75
Preservative	None
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C or -80°C. Avoid repeated freeze / thaw cycles. Preservative: None Constituents: 0.001% Tween 20, 30mM HEPES, 2mM EDTA, 150mM Sodium chloride, pH 6.75

GENE INFORMATION

Gene Name	CAPN11 calpain 11 [Homo sapiens (human)]
Official Symbol	CAPN11
Synonyms	CAPN11; calpain 11; calpain11; calpain-11; CANP 11; calcium-activated neutral proteinase 11;

Entrez Gene ID	11131
mRNA Refseq	NM_007058.3
Protein Refseq	NP_008989.2
UniProt ID	Q9UMQ6
Chromosome Location	6p12
Pathway	Integrin-mediated cell adhesion, organism-specific biosystem;
Function	calcium ion binding; calcium-dependent cysteine-type endopeptidase activity; peptidase activity;