



Human CCHCR1 peptide (DAG-P0318)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Antigen Description	CCHCR1 (Coiled coil alpha helical rod protein 1) is a regulator of keratinocyte proliferation or differentiation. It is overexpressed in keratinocytes in psoriatic lesions.
Purity	70 - 90% by HPLC.
Conjugate	Unconjugated
Format	Liquid
Preservative	None
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C or -80°C. Avoid repeated freeze / thaw cycles. Information available upon request.

GENE INFORMATION

Gene Name	CCHCR1 coiled-coil alpha-helical rod protein 1 [Homo sapiens (human)]
Official Symbol	CCHCR1
Synonyms	CCHCR1; coiled-coil alpha-helical rod protein 1; HCR; SBP; C6orf18; pg8; StAR-binding protein; putative gene 8 protein; alpha-helical coiled-coil rod protein; HCR (a-helix coiled-coil rod homologue);
Entrez Gene ID	54535
mRNA Refseq	NM_001105563.1
Protein Refseq	NP_001099033.1
UniProt ID	Q2TB68

