



Human DPPA4 blocking peptide (DAG-P0436)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Antigen Description	May be involved in the maintenance of active epigenetic status of target genes (By similarity). May inhibit differentiation of embryonic cells into a primitive ectoderm lineage.
Conjugate	Unconjugated
Applications	BL
Format	Liquid
Preservative	None
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C or -80°C. Avoid repeated freeze / thaw cycles. Information available upon request.

GENE INFORMATION

Gene Name	DPPA4 developmental pluripotency associated 4 [Homo sapiens (human)]
Official Symbol	DPPA4
Synonyms	DPPA4; developmental pluripotency associated 4; 2410091M23Rik; developmental pluripotency-associated protein 4;
Entrez Gene ID	55211
mRNA Refseq	NM_018189.3
Protein Refseq	NP_060659.3
UniProt ID	Q7L190
Chromosome Location	3q13.13

Function	protein binding;
----------	------------------
