



Human ETV2 peptide (DAG-P0492)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Antigen Description	ETV2 is a member of the ETS family of transcription factors. It has been proposed that ETV2 acts downstream of BMP, notch and wnt signalling in blood and vessel specification. ER71 has also shown to be important for Flk1 expression as well as hematopoietic and vascular development.
Purity	70 - 90% by HPLC.
Conjugate	Unconjugated
Format	Liquid
Preservative	None
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C or -80°C. Avoid repeated freeze / thaw cycles. Information available upon request.

GENE INFORMATION

Gene Name	ETV2 ets variant 2 [Homo sapiens (human)]
Official Symbol	ETV2
Synonyms	ETV2; ets variant 2; ER71; ETSRP71; ETS translocation variant 2; ets variant gene 2; ets-related protein 71;
Entrez Gene ID	2116
mRNA Refseq	NM_014209.2
Protein Refseq	NP_055024.2
UniProt ID	B3KUL0

Chromosome Location	19q13.12
Function	sequence-specific DNA binding; sequence-specific DNA binding RNA polymerase II transcription factor activity;