



Human C19ORF40 peptide (DAG-P0499)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Antigen Description	FAAP24 is a component of the Fanconi anemia (FA) core complex (see MIM 227650), which plays a crucial role in DNA damage response (Ciccia et al., 2007 [PubMed 17289582]).[supplied by OMIM, Mar 2008]
Conjugate	Unconjugated
Format	Liquid
Preservative	None
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C or -80°C. Avoid repeated freeze / thaw cycles. Information available upon request.

GENE INFORMATION

Gene Name	C19orf40 chromosome 19 open reading frame 40 [Homo sapiens (human)]
Official Symbol	C19ORF40
Synonyms	C19ORF40; chromosome 19 open reading frame 40; FAAP24; Fanconi anemia-associated protein of 24 kDa; Fanconi anemia-associated protein, 24 kDa;
Entrez Gene ID	91442
mRNA Refseq	NM_152266.3
Protein Refseq	NP_689479.1
UniProt ID	Q9BTP7
Chromosome Location	19q13.11

Pathway	DNA Repair, organism-specific biosystem; FA core complex, organism-specific biosystem; FA core complex, conserved biosystem; Fanconi Anemia pathway, organism-specific biosystem; Fanconi anemia pathway, organism-specific biosystem; Fanconi anemia pathway, conserved biosystem;
Function	DNA binding; chromatin binding; protein binding;