



Human HNRNPA2B1 peptide (DAG-P0587)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Antigen Description	This gene belongs to the A/B subfamily of ubiquitously expressed heterogeneous nuclear ribonucleoproteins (hnRNPs). The hnRNPs are RNA binding proteins and they complex with heterogeneous nuclear RNA (hnRNA). These proteins are associated with pre-mRNAs in the nucleus and appear to influence pre-mRNA processing and other aspects of mRNA metabolism and transport. While all of the hnRNPs are present in the nucleus, some seem to shuttle between the nucleus and the cytoplasm. The hnRNP proteins have distinct nucleic acid binding properties. The protein encoded by this gene has two repeats of quasi-RRM domains that bind to RNAs. This gene has been described to generate two alternatively spliced transcript variants which encode different isoforms. [provided by RefSeq, Jul 2008]
Conjugate	Unconjugated
Sequence Similarities	Contains 2 RRM (RNA recognition motif) domains.
Format	Lyophilized powder
Preservative	None
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C or -80°C. Avoid repeated freeze / thaw cycles. Information available upon request.

GENE INFORMATION

Gene Name	HNRNPA2B1 heterogeneous nuclear ribonucleoprotein A2/B1 [Homo sapiens (human)]
Official Symbol	HNRNPA2B1
Synonyms	HNRNPA2B1; heterogeneous nuclear ribonucleoprotein A2/B1; RNPA2; HNRPA2; HNRPB1; SNRPB1; HNRNPA2; HNRNPB1; IBMPFD2; HNRPA2B1; heterogeneous nuclear ribonucleoproteins A2/B1; hnRNP A2 / hnRNP B1; nuclear ribonucleoprotein particle A2 protein;

Entrez Gene ID	3181
mRNA Refseq	NM_002137.3
Protein Refseq	NP_002128.1
UniProt ID	P22626
Chromosome Location	7p15
Pathway	Gene Expression, organism-specific biosystem; Processing of Capped Intron-Containing Pre-mRNA, organism-specific biosystem; mRNA Splicing, organism-specific biosystem; mRNA Splicing - Major Pathway, organism-specific biosystem; mRNA processing, organism-specific biosystem;
Function	RNA binding; nucleotide binding; poly(A) RNA binding; pre-mRNA intronic binding; protein binding; single-stranded telomeric DNA binding;