



Human JMJD4 peptide (DAG-P0604)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Conjugate	Unconjugated
Sequence Similarities	Contains 1 JmjC domain.
Format	Liquid
Preservative	None
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C or -80°C. Avoid repeated freeze / thaw cycles. Information available upon request.

GENE INFORMATION

Gene Name	JMJD4 jumonji domain containing 4 [Homo sapiens (human)]
Official Symbol	JMJD4
Synonyms	JMJD4; jumonji domain containing 4; jmjC domain-containing protein 4; jumonji domain-containing protein 4;
Entrez Gene ID	65094
mRNA Refseq	NM_001161465.1
Protein Refseq	NP_001154937.1
UniProt ID	Q9H9V9
Chromosome Location	1q42.13
Function	protein binding;