



Human HPC2 peptide (DAG-P0649)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Antigen Description	HPC2 is a member of the HIR complex that also contains HIR1, HIR2 and HIR3. The HIR complex is responsible for the repression of histone gene transcription as well having a role in nucleosome assembly.
Purity	70 - 90% by HPLC.
Conjugate	Unconjugated
Format	Liquid
Preservative	None
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C or -80°C. Avoid repeated freeze / thaw cycles. Information available upon request.

GENE INFORMATION

Gene Name	HPC2 Hpc2p [Saccharomyces cerevisiae S288c]
Synonyms	Hpc2p;
Entrez Gene ID	852516
mRNA Refseq	NM_001178563.3
Protein Refseq	NP_009774.4
UniProt ID	Q01448
Chromosome Location	chromosome: II
Function	contributes_to DNA binding; contributes_to nucleosome binding;