



MKI67 blocking peptide (DAG-P0692)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Antigen Description	This gene encodes a nuclear protein that is associated with and may be necessary for cellular proliferation. Alternatively spliced transcript variants have been described. A related pseudogene exists on chromosome X. [provided by RefSeq, Mar 2009]
Conjugate	Unconjugated
Applications	BL
Sequence Similarities	Contains 1 FHA domain.
Format	Liquid
Preservative	None
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C or -80°C. Avoid repeated freeze / thaw cycles. Information available upon request.

GENE INFORMATION

Gene Name	MKI67 marker of proliferation Ki-67 [Homo sapiens (human)]
Official Symbol	MKI67
Synonyms	MKI67; marker of proliferation Ki-67; KIA; MIB-1; antigen Ki-67; proliferation-related Ki-67 antigen; antigen identified by monoclonal antibody Ki-67;
Entrez Gene ID	4288
mRNA Refseq	NM_001145966.1
Protein Refseq	NP_001139438.1

UniProt ID	P46013
Chromosome Location	10q26.2
Function	ATP binding; poly(A) RNA binding; protein C-terminus binding; protein binding;