



Mouse LHX8 peptide (DAG-P0734)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Antigen Description	Transcription factor involved in differentiation of certain neurons and mesenchymal cells.
Conjugate	Unconjugated
Sequence Similarities	Contains 1 homeobox DNA-binding domain.Contains 2 LIM zinc-binding domains.
Format	Liquid
Preservative	None
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C or -80°C. Avoid repeated freeze / thaw cycles. Information available upon request.

GENE INFORMATION

Gene Name	Lhx8 LIM homeobox protein 8 [Mus musculus (house mouse)]
Official Symbol	LHX8
Synonyms	LHX8; LIM homeobox protein 8; L3; Lhx7; LIM/homeobox protein Lhx8; LIM homeobox protein 7; LIM homeo box protein 8; LIM/homeobox protein Lhx7;
Entrez Gene ID	16875
mRNA Refseq	NM_010713.2
Protein Refseq	NP_034843.2
UniProt ID	O35652
Chromosome Location	3; 3 H3-H4

Function

DNA binding; metal ion binding; sequence-specific DNA binding; sequence-specific DNA binding transcription factor activity; zinc ion binding;
