



Human LSM11 peptide (DAG-P0771)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Conjugate	Unconjugated
Format	Liquid
Preservative	None
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C or -80°C. Avoid repeated freeze / thaw cycles. Information available upon request.

GENE INFORMATION

Gene Name	LSM11 LSM11, U7 small nuclear RNA associated [Homo sapiens (human)]
Official Symbol	LSM11
Synonyms	LSM11; LSM11, U7 small nuclear RNA associated; U7 snRNA-associated Sm-like protein LSm11; 2210404M20Rik; LSM11 homolog, U7 small nuclear RNA associated;
Entrez Gene ID	134353
mRNA Refseq	NM_173491.2
Protein Refseq	NP_775762.1
UniProt ID	P83369
Chromosome Location	5q33.3
Pathway	Cleavage of Growing Transcript in the Termination Region, organism-specific biosystem; Gene Expression, organism-specific biosystem; Processing of Capped Intronless Pre-mRNA, organism-specific biosystem; RNA Polymerase II Transcription, organism-specific biosystem; RNA Polymerase II Transcription Termination, organism-specific biosystem; SLBP Dependent

Processing of Replication-Dependent Histone Pre-mRNAs, organism-specific biosystem; SLBP independent Processing of Histone Pre-mRNAs, organism-sp

Function	U7 snRNA binding; protein binding;
-----------------	------------------------------------
