



Human LAMP2 blocking peptide (DAG-P0775)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Antigen Description	The protein encoded by this gene is a member of a family of membrane glycoproteins. This glycoprotein provides selectins with carbohydrate ligands. It may play a role in tumor cell metastasis. It may also function in the protection, maintenance, and adhesion of the lysosome. Alternative splicing of this gene results in multiple transcript variants encoding distinct proteins. [provided by RefSeq, Jul 2008]
Conjugate	Unconjugated
Applications	BL
Format	Liquid
Preservative	None
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C or -80°C. Avoid repeated freeze / thaw cycles. Information available upon request.

GENE INFORMATION

Gene Name	LAMP2 lysosomal-associated membrane protein 2 [Homo sapiens (human)]
Official Symbol	LAMP2
Synonyms	LAMP2; lysosomal-associated membrane protein 2; LAMPB; CD107b; LAMP-2; LGP110; lysosome-associated membrane glycoprotein 2; CD107 antigen-like family member B;
Entrez Gene ID	3920
mRNA Refseq	NM_001122606.1
Protein Refseq	NP_001116078.1

UniProt ID	P13473
Chromosome Location	Xq24
Pathway	Hemostasis, organism-specific biosystem; Lysosome, organism-specific biosystem; Lysosome, conserved biosystem; Phagosome, organism-specific biosystem; Phagosome, conserved biosystem; Platelet activation, signaling and aggregation, organism-specific biosystem; Platelet degranulation, organism-specific biosystem; Response to elevated platelet cytosolic Ca ²⁺ , organism-specific biosystem; Senescence and Autophagy, organism-specific biosystem; Tuberculosis, organism-specific biosystem; Tuberculosis,
Function	enzyme binding;