



MMP24 peptide (DAG-P0802)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Antigen Description	Proteins of the matrix metalloproteinase (MMP) family are involved in the breakdown of extracellular matrix in normal physiological processes, such as embryonic development, reproduction, and tissue remodeling, as well as in disease processes, such as arthritis and metastasis. Most MMPs are secreted as inactive proproteins which are activated when cleaved by extracellular proteinases. However, the protein encoded by this gene is a member of the membrane-type MMP (MT-MMP) subfamily; each member of this subfamily contains a potential transmembrane domain suggesting that these proteins are expressed at the cell surface rather than secreted. This protein activates MMP2 by cleavage. The gene has previously been referred to as MMP25 but has been renamed MMP24. [provided by RefSeq, Jul 2008]
Specificity	Predominantly expressed in brain, kidney, pancreas and lung. Overexpressed in a series of brain tumors, including astrocytomas and glioblastomas.
Purity	> 95 % by SDS-PAGE.
Conjugate	Unconjugated
Applications	WB, ELISA
Sequence Similarities	Belongs to the peptidase M10A family. Contains 4 hemopexin-like domains.
Format	Liquid
Buffer	Preservative: None Constituents: 0.001% Tween 20, 30mM HEPES, 2mM EDTA, 150mM Sodium chloride, pH 6.75
Preservative	None
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C or -80°C. Avoid repeated freeze / thaw cycles. Preservative: None Constituents: 0.001% Tween 20, 30mM HEPES, 2mM EDTA, 150mM Sodium chloride, pH 6.75

GENE INFORMATION

Gene Name	MMP24 matrix metalloproteinase 24 (membrane-inserted) [Homo sapiens (human)]
Official Symbol	MMP24
Synonyms	MMP24; matrix metalloproteinase 24 (membrane-inserted); MMP25; MMP-24; MT5MMP; MTMMP5; MT-MMP5; MT5-MMP; MT-MMP 5; matrix metalloproteinase-24; membrane-type 5 matrix metalloproteinase; membrane-type matrix metalloproteinase 5; membrane-type-5 matrix metalloproteinase; matrix metalloproteinase 24 (membrane-inserted);
Entrez Gene ID	10893
mRNA Refseq	NM_006690.3
Protein Refseq	NP_006681.1
UniProt ID	Q86VV6
Chromosome Location	20q11.2
Pathway	Activation of Matrix Metalloproteinases, organism-specific biosystem; Degradation of the extracellular matrix, organism-specific biosystem; Extracellular matrix organization, organism-specific biosystem; Matrix Metalloproteinases, organism-specific biosystem;
Function	calcium ion binding; enzyme activator activity; metalloendopeptidase activity; zinc ion binding;