



Mouse NTNG2 peptide (DAG-P0969)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Antigen Description	Promotes neurite outgrowth of both axons and dendrites.
Conjugate	Unconjugated
Sequence Similarities	Contains 3 laminin EGF-like domains.Contains 1 laminin N-terminal domain.
Format	Liquid
Preservative	None
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C or -80°C. Avoid repeated freeze / thaw cycles. Information available upon request.

GENE INFORMATION

Gene Name	NTNG2 netrin G2 [Homo sapiens (human)]
Official Symbol	NTNG2
Synonyms	NTNG2; netrin G2; Lmnt2; NTNG1; LHLL9381; bA479K20.1; netrin-G2; laminet 2; laminet-2; netrin G1; bA479K20.1 (novel protein);
Entrez Gene ID	84628
mRNA Refseq	NM_032536.2
Protein Refseq	NP_115925.2
UniProt ID	O95898
Chromosome Location	9q34

Pathway	Axon guidance, organism-specific biosystem; Axon guidance, conserved biosystem; Cell adhesion molecules (CAMs), organism-specific biosystem; Cell adhesion molecules (CAMs), conserved biosystem;
Function	molecular_function; protein binding;