



Human RAP2A peptide (DAG-P1110)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Antigen Description	Rap2 is a 21 kDa GTP-binding protein related to Ras. It is implicated in vesiculation and synaptic plasticity.
Conjugate	Unconjugated
Format	Liquid
Preservative	None
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C or -80°C. Avoid repeated freeze / thaw cycles. Information available upon request.

GENE INFORMATION

Gene Name	RAP2A RAP2A, member of RAS oncogene family [Homo sapiens (human)]
Official Symbol	RAP2A
Synonyms	RAP2A; RAP2A, member of RAS oncogene family; KREV; RAP2; K-REV; RbBP-30; ras-related protein Rap-2a; RAP2, member of RAS oncogene family (K-rev);
Entrez Gene ID	5911
mRNA Refseq	NM_021033.6
Protein Refseq	NP_066361.1
UniProt ID	P10114
Chromosome Location	13q34
Pathway	B Cell Receptor Signaling Pathway, organism-specific biosystem;

Function

GTP binding; GTPase activity; protein binding;
