



Human SIRT3 peptide (DAG-P1120)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Antigen Description	This gene encodes a member of the sirtuin family of proteins, homologs to the yeast Sir2 protein. Members of the sirtuin family are characterized by a sirtuin core domain and grouped into four classes. The functions of human sirtuins have not yet been determined; however, yeast sirtuin proteins are known to regulate epigenetic gene silencing and suppress recombination of rDNA. Studies suggest that the human sirtuins may function as intracellular regulatory proteins with mono-ADP-ribosyltransferase activity. The protein encoded by this gene is included in class I of the sirtuin family. Two alternatively spliced transcript variants that encode different proteins have been described for this gene. [provided by RefSeq, Jul 2008]
Specificity	Widely expressed.
Conjugate	Unconjugated
Sequence Similarities	Belongs to the sirtuin family. Class I subfamily. Contains 1 deacetylase sirtuin-type domain.
Preservative	None
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C or -80°C. Avoid repeated freeze / thaw cycles. Information available upon request.

GENE INFORMATION

Gene Name	SIRT3 sirtuin 3 [Homo sapiens (human)]
Official Symbol	SIRT3
Synonyms	SIRT3; sirtuin 3; SIR2L3; NAD-dependent protein deacetylase sirtuin-3, mitochondrial; sir2-like 3; sirtuin type 3; SIR2-like protein 3; regulatory protein SIR2 homolog 3; NAD-dependent deacetylase sirtuin-3, mitochondrial; silent mating type information regulation 2, S.cerevisiae, homolog 3; mitochondrial nicotinamide adenine dinucleotide-dependent deacetylase;

Entrez Gene ID	23410
mRNA Refseq	NM_001017524.2
Protein Refseq	NP_001017524.1
UniProt ID	Q9NTG7
Chromosome Location	11p15.5
Pathway	Energy Metabolism, organism-specific biosystem; Signaling events mediated by HDAC Class I, organism-specific biosystem; Signaling events mediated by HDAC Class III, organism-specific biosystem;
Function	NOT NAD+ ADP-ribosyltransferase activity; NAD+ binding; NAD-dependent histone deacetylase activity (H3-K14 specific); protein binding; zinc ion binding;