



Human ARHGAP24 peptide (DAG-P1129)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Antigen Description	ARHGAPs, such as ARHGAP24, encode negative regulators of Rho GTPases (see ARHA; MIM 165390), which are implicated in actin remodeling, cell polarity, and cell migration (Katoh and Katoh, 2004 [PubMed 15254788]).[supplied by OMIM, Mar 2008]
Purity	70 - 90% by HPLC.
Conjugate	Unconjugated
Format	Liquid
Preservative	None
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C or -80°C. Avoid repeated freeze / thaw cycles. Information available upon request.

GENE INFORMATION

Gene Name	ARHGAP24 Rho GTPase activating protein 24 [Homo sapiens (human)]
Official Symbol	ARHGAP24
Synonyms	ARHGAP24; Rho GTPase activating protein 24; p73; FILGAP; RCGAP72; RC-GAP72; p73RhoGAP; rho GTPase-activating protein 24; rhoGAP of 73 kDa; sarcoma antigen NY-SAR-88; filamin-A-associated RhoGAP; rho-type GTPase-activating protein 24; RAC1- and CDC42-specific GTPase-activating protein of 72 kDa;
Entrez Gene ID	83478
mRNA Refseq	NM_001025616.2
Protein Refseq	NP_001020787.2

UniProt ID	Q8N264
Chromosome Location	4q22.1
Pathway	Rho GTPase cycle, organism-specific biosystem; Signal Transduction, organism-specific biosystem; Signaling by Rho GTPases, organism-specific biosystem;
Function	Rac GTPase activator activity; protein binding;