



Human SMCHD1 peptide (DAG-P1149)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Antigen Description	This gene encodes a protein which contains a hinge region domain found in members of the SMC (structural maintenance of chromosomes) protein family. [provided by RefSeq, Dec 2011]
Conjugate	Unconjugated
Format	Liquid
Preservative	None
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C or -80°C. Avoid repeated freeze / thaw cycles. Information available upon request.

GENE INFORMATION

Gene Name	SMCHD1 structural maintenance of chromosomes flexible hinge domain containing 1 [Homo sapiens (human)]
Official Symbol	SMCHD1
Synonyms	SMCHD1; structural maintenance of chromosomes flexible hinge domain containing 1; structural maintenance of chromosomes flexible hinge domain-containing protein 1; SMC hinge domain-containing protein 1;
Entrez Gene ID	23347
mRNA Refseq	NM_015295.2
Protein Refseq	NP_056110.2
UniProt ID	A6NHR9
Chromosome Location	18p11.32

Function

ATP binding;