



## Mouse STRA8 peptide (DAG-P1202)

This product is for research use only and is not intended for diagnostic use.

### PRODUCT INFORMATION

<b>Antigen Description</b>	Required for the transition into meiosis for both female and male germ cells. In female germ cells, required for premeiotic DNA replication and subsequent events in meiotic prophase.
<b>Specificity</b>	Expressed specifically in testis.
<b>Conjugate</b>	Unconjugated
<b>Format</b>	Liquid
<b>Preservative</b>	None
<b>Storage</b>	Shipped at 4°C. Upon delivery aliquot and store at -20°C or -80°C. Avoid repeated freeze / thaw cycles. Information available upon request.

### GENE INFORMATION

<b>Gene Name</b>	<a href="#">Stra8 stimulated by retinoic acid gene 8 [ Mus musculus (house mouse) ]</a>
<b>Official Symbol</b>	STRA8
<b>Synonyms</b>	STRA8; stimulated; stimulated by retinoic acid gene 8 protein;
<b>Entrez Gene ID</b>	<a href="#">20899</a>
<b>mRNA Refseq</b>	<a href="#">NM_009292.1</a>
<b>Protein Refseq</b>	<a href="#">NP_033318.1</a>
<b>UniProt ID</b>	P70278
<b>Chromosome Location</b>	6 B1; 6

Function	protein binding;
----------	------------------

---