



Mouse STRA8 peptide (DAG-P1203)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Antigen Description	Required for the transition into meiosis for both female and male germ cells. In female germ cells, required for premeiotic DNA replication and subsequent events in meiotic prophase.
Specificity	Expressed specifically in testis.
Conjugate	Unconjugated
Format	Liquid
Preservative	None
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C or -80°C. Avoid repeated freeze / thaw cycles. Information available upon request.

GENE INFORMATION

Gene Name	Stra8 stimulated by retinoic acid gene 8 [Mus musculus (house mouse)]
Official Symbol	STRA8
Synonyms	STRA8; stimulated; stimulated by retinoic acid gene 8 protein;
Entrez Gene ID	20899
mRNA Refseq	NM_009292.1
Protein Refseq	NP_033318.1
UniProt ID	P70278
Chromosome Location	6 B1; 6

Function

protein binding;