



Human ZSCAN10 peptide (DAG-P1292)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Antigen Description	Znf206 is a zinc finger protein and probable transcriptional regulator. It is expressed in embryonic stem cells and downregulated upon differentiation. Both Oct4 and Sox2 bind the znf206 gene.
Conjugate	Unconjugated
Format	Liquid
Preservative	None
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C or -80°C. Avoid repeated freeze / thaw cycles. Information available upon request.

GENE INFORMATION

Gene Name	ZSCAN10 zinc finger and SCAN domain containing 10 [Homo sapiens (human)]
Official Symbol	ZSCAN10
Synonyms	ZSCAN10; zinc finger and SCAN domain containing 10; ZNF206; zinc finger and SCAN domain-containing protein 10; zinc finger protein 206;
Entrez Gene ID	84891
mRNA Refseq	NM_001282415.1
Protein Refseq	NP_001269344.1
UniProt ID	Q96SZ4
Chromosome Location	16p13.3

Function	metal ion binding; sequence-specific DNA binding; sequence-specific DNA binding transcription factor activity; sequence-specific DNA binding transcription factor activity;
-----------------	---
