



Human ZNF217 peptide (DAG-P1293)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Antigen Description	Binds to the promoters of target genes and functions as repressor. Promotes cell proliferation and antagonizes cell death. Promotes phosphorylation of AKT1 at 'Ser-473'.
Conjugate	Unconjugated
Sequence Similarities	Belongs to the krueppel C2H2-type zinc-finger protein family. Contains 7 C2H2-type zinc fingers.
Format	Liquid
Preservative	None
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C or -80°C. Avoid repeated freeze / thaw cycles. Information available upon request.

GENE INFORMATION

Gene Name	ZNF217 zinc finger protein 217 [Homo sapiens (human)]
Official Symbol	ZNF217
Synonyms	ZNF217; zinc finger protein 217; ZABC1;
Entrez Gene ID	7764
mRNA Refseq	NM_006526.2
Protein Refseq	NP_006517.1
UniProt ID	O75362
Chromosome Location	20q13.2

Function

metal ion binding; protein binding; sequence-specific DNA binding transcription factor activity;
transcription regulatory region DNA binding;
