



Human RCOR2 peptide (DAG-P1443)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Antigen Description	May act as a component of a corepressor complex that represses transcription.
Conjugate	Unconjugated
Sequence Similarities	Belongs to the CoREST family. Contains 1 ELM2 domain. Contains 2 SANT domains.
Format	Liquid
Preservative	None
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C or -80°C. Avoid repeated freeze / thaw cycles. Information available upon request.

GENE INFORMATION

Gene Name	RCOR2 REST corepressor 2 [Homo sapiens (human)]
Official Symbol	RCOR2
Synonyms	RCOR2; REST corepressor 2;
Entrez Gene ID	283248
mRNA Refseq	NM_173587.3
Protein Refseq	NP_775858.2
UniProt ID	Q8IZ40
Chromosome Location	11q13.1
Function	DNA binding; chromatin binding; sequence-specific DNA binding transcription factor activity;

transcription corepressor activity;
