



Human NAA25 peptide (DAG-P1779)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Antigen Description	This gene encodes the auxiliary subunit of the heteromeric N-terminal acetyltransferase B complex. This complex acetylates methionine residues that are followed by acidic or asparagine residues.[provided by RefSeq, Mar 2010]
Purity	70 - 90% by HPLC.
Conjugate	Unconjugated
Sequence Similarities	Belongs to the MDM20/NAA25 family.Contains 4 TPR repeats.
Format	Liquid
Preservative	None
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C or -80°C. Avoid repeated freeze / thaw cycles. Information available upon request.

GENE INFORMATION

Gene Name	NAA25 N(alpha)-acetyltransferase 25, NatB auxiliary subunit [Homo sapiens (human)]
Official Symbol	NAA25
Synonyms	NAA25; N(alpha)-acetyltransferase 25, NatB auxiliary subunit; NAP1; MDM20; C12orf30; N-alpha-acetyltransferase 25, NatB auxiliary subunit; natB complex subunit MDM20; mitochondrial distribution and morphology 20; N-terminal acetyltransferase B complex subunit NAA25;
Entrez Gene ID	80018
mRNA Refseq	NM_024953.3

Protein Refseq	NP_079229.2
UniProt ID	Q14CX7
Chromosome Location	12q24.13