



PGRMC2 peptide (DAG-P1832)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Antigen Description	Membrane associated progesterone receptor component 2 (PGRMC2) belongs to the cytochrome b5 family. The molecular function is primarily steroid binding.
Conjugate	Unconjugated
Format	Lyophilised
Preservative	None
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid freeze / thaw cycles.

GENE INFORMATION

Gene Name	PGRMC2 progesterone receptor membrane component 2 [Homo sapiens (human)]
Official Symbol	PGRMC2
Synonyms	PGRMC2; progesterone receptor membrane component 2; DG6; PMBP; membrane-associated progesterone receptor component 2; steroid receptor protein DG6; progesterone membrane binding protein;
Entrez Gene ID	10424
mRNA Refseq	NM_006320.4
Protein Refseq	NP_006311.2
UniProt ID	O15173
Chromosome Location	4q26
Function	heme binding; steroid binding; steroid hormone receptor activity;