



Prickle peptide (DAG-P1842)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Conjugate	Unconjugated
Format	Liquid
Preservative	None
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C or -80°C. Avoid repeated freeze / thaw cycles. Information available upon request.

BACKGROUND

Introduction	Prickle is a LIM domain protein that appears to serve as a nuclear receptor for REST/NRSF, REST4, and possibly other transcription factors. The human protein is a 831 amino acid protein containing 3 N-terminal LIM domains and 3 C-terminal nuclear localiza
Keywords	FLJ31627; FLJ31937; MGC138902; MGC138903; Prickle 1; Prickle homolog 1 (Drosophila); Prickle like protein 1; Prickle-like 1; Prickle1; REST-Interacting LIM Domain Protein; REST/NRSF interacting LIM domain protein 1; RILP; Prickle